

equinox[®]
PLATFORM SHOULDER SYSTEM




EQUINOXE[®] REVERSE

ANATOMICAL. REDEFINED.

**7x Reduction
 in Scapular Notching¹
 ...and no sacrifice of glenoid fixation or stability.²**

**Reported Scapular Notching Rate for
 Grammont-Style Reverse Shoulder Prosthesis**




Grammont-Style Prosthesis, weighted average of 8 studies, n=868¹
 Equinox Reverse, 7 site multicenter study, n=226¹



Minimized Scapular Notching

Lateralized Humeral Components

- The three larger diameter glenospheres (38, 42 and 46mm) lateralize the humerus (without lateralizing the center of rotation) and increase joint stability.^{3,4}



**Grammont
Style**

Equinoxe


- Decreasing the humeral neck angle to 145 degrees further lateralizes the humerus (without lateralizing the center of rotation) and helps restore anatomic tensioning of the remaining intact rotator cuff.^{3,4} Additionally, placing the humeral tray on top of the resection eliminates the need to conically ream the proximal humerus, improves exposure and allows for larger glenospheres to be implanted (i.e., the size of the proximal humerus does not dictate the size of the glenosphere).

Greater Range of Motion

- The innovative glenoid baseplate design has a built-in offset which distally shifts the glenosphere to a position that prevents humeral liner impingement on the inferior glenoid. This offset negates the need for additional bone-consuming implantation techniques (i.e., inferiorly tilting the baseplate or pre-notching the bone).^{3,4}
- The increased stability provided by the larger diameter glenospheres enable the humeral liners to be less constrained relative to other systems and thereby permits greater range of motion prior to impingement.^{3,4}
- The extended glenosphere articular surface and chamfered sides maximize inferior overhang designed to minimize the potential for scapular notching and improve range of motion.




**Innovative
Baseplate**



**Larger
Glenosphere**

equinox[®]
PLATFORM SHOULDER SYSTEM



**Structural Integrity.
Even in the Most
Challenging Conditions.**



Enhanced Glenoid Fixation

Minimized Torque on the Glenoid

- The larger diameter glenospheres result in a medialized center of rotation, thereby minimizing the torque on the glenoid. This medialized center of rotation has a long history of clinical success.⁵
- The curved back of the glenoid baseplate also allows for fixation by converting destabilizing shear forces into stabilizing compressive forces.



Stable Construct

Strong Initial Fixation

- Strong initial fixation can be achieved with the press-fit bone cage of the glenoid baseplate, while the six-hole baseplate design provides up to 30 degrees of screw variability to ensure optimal compression screw placement and purchase, even in poor quality bone.⁶
- Locking caps are provided to secure the compression screws to the glenoid baseplate at the desired variable angle.
- Bench testing conducted on the Equinox[®] reverse quantified micromotion values at approximately half of those published with other systems measuring fixation using similar testing methodologies.⁶⁻⁹




Multiple Options for Screw Placement


Long-Term Biologic Fixation

- Unique to the Equinox, bone graft can be inserted into the cage to promote bone through-growth,¹⁰ which enhances the probability of long-term biologic fixation.

Allows for Bone Through-Growth¹⁰



Unique Bone Cage



equinox[®]
PLATFORM SHOULDER SYSTEM

One Stem.
Two Options.



Seamlessly Integrated System




**Revision of a TSA
to a Reverse**

Standardized Humeral Preparation

- The Equinox's platform system enables a surgeon to convert from a total shoulder to a reverse without humeral stem removal.
- Using the same humeral stem, humeral instrumentation and humeral osteotomy for both indications standardizes the procedure and empowers the surgeon to intra-operatively decide 'primary vs. reverse'.

Standardized Glenoid Preparation

- The offset bone cage of the glenoid baseplate is sized and positioned to be placed in the center of the glenoid to fill a central bone defect while distally shifting the glenosphere to ensure inferior overhang.
- The six screw holes of the glenoid baseplate are positioned to provide screw fixation, even when revising a pegged or keeled glenoid to a reverse shoulder.



**Revision
Potential of
Pegged and
Keeled Glenoid**

Conclusion

Scapular notching is currently addressed in the marketplace with implant designs that either a) lateralize the center of rotation, which causes greater torque on the glenoid, or b) require additional bone-consuming surgical techniques such as inferiorly tilting the baseplate or pre-notching the bone.¹¹ Both of these options create the potential for long-term glenoid fixation

challenges. The Equinox Reverse Shoulder, however, minimizes scapular notching exclusively in design while maintaining a medialized center of rotation without bone consuming techniques. These critical attributes, along with the seamlessly integrated platform stem, differentiate the Equinox and provide a compelling reason for a surgeon to *Experience the Power of the Equinox.*

References

1. **Roche C, Wright T, Flurin PH, Grey S, Jones R, Routman H, Gilot G, Zuckerman JD.** Scapular notching radiographic analysis: recommendations for glenoid plate positioning and glenosphere overhang in male and female patients. Transactions of the 58th Annual Meeting of the Orthopaedic Research Society; 2012 Feb 4-7; San Francisco, CA.
2. **Flurin P, et al.** A correlation of five commonly used clinical metrics to measure outcomes in shoulder arthroplasty. Poster presentation at the 58th Annual Meeting of the Orthopaedic Research Society; 2012 Feb 4-7; San Francisco, CA.
3. **Roche C, Flurin PH, Wright T, Crosby LA, Mauldin M, Zuckerman JD.** Geometric analysis of the grammont reverse shoulder prosthesis: an evaluation of the relationship between prosthetic design parameters and clinical failure modes. Proceedings of the 19th Annual Congress of the International Society for Technology in Arthroplasty; 2006 Oct 6-9; New York, NY.
4. **Roche C, Flurin PH, Wright T, Crosby LA, Mauldin M, Zuckerman JD.** An evaluation of the relationships between reverse shoulder design parameters and range of motion, impingement, and stability. *J Shoulder Elbow Surg.* 2009 Sept-Oct;18(5):734-41.
5. **Boileau P, Watkinson DJ, Hatzidakis AM, Balg F.** Grammont reverse prosthesis: design, rationale and biomechanics. *J Shoulder Elbow Surg.* 2005 Jan-Feb;14(1 Suppl S):147S-161S.
6. **Roche C, Flurin PH, Wright T, Crosby LA, Zuckerman JD.** Effect of varying angle screw configuration and bone density on reverse shoulder glenoid fixation following cyclic loading. Transactions of the 54th Annual Orthopaedic Research Society Meeting; 2008 Mar 2-5; San Francisco, CA.
7. **Harman M, Frankle M, Vasey M, Banks S.** Initial glenoid component fixation in "reverse" total shoulder arthroplasty: a biomechanical evaluation. *J Shoulder Elbow Surg.* 2005 Jan-Feb;14(1 Suppl S):162S-167S.
8. **Virani N, et al.** Experimental and finite element analysis of baseplate micromotion in the reverse shoulder design. Transactions of the 53rd Annual Meeting of the Orthopaedic Research Society; 2007 Feb 11-14; San Diego, CA.
9. Data on file at Exactech.
10. Animal study data on file at Exactech.
11. **Nyffeler RW, Werner CM, Gerber C.** Biomechanical relevance of glenoid component positioning in the reverse Delta III total shoulder prosthesis. *J Shoulder Elbow Surg.* 2005 Sep-Oct;14(5):524-8.

352-377-1140
1-800-EXACTECH
www.exac.com

Exactech®

718-04-90 Rev. A
Equinox Reverse Campaign Brochure 0712

