



InterSep™
Calcium Sulfate

Surgeon focused. Patient driven.™ **Exactech®**

InterSep™

Calcium Sulfate

Continuing our goal to provide infection-related solutions, Exactech introduces InterSep – a 100 percent synthetic calcium sulfate bone void filler engineered to fully resorb and replace bone during the healing process. InterSep may be used in an infected bone site and provides surgeon flexibility with its bead and paste options.

ADVANTAGES

- May be used in an infected site
- Provides case-by-case flexibility by offering various bead sizes (3mm, 4.5mm and 6mm) and a paste option
- Bone graft substitute that resorbs and is replaced with bone during the healing process
- No impurities, pure synthetic material
- Biodegradable and biocompatible
- Easy to mix and deliver
- Fast setting time
- Radiopaque



Product	Product Codes	Paste Volume	Bead Volume*			Setting Time
			3mm	4.5mm	6mm	
InterSep™ Calcium Sulfate	CS-05CC	5cc	9cc	11cc	12.5cc	5 min
	CS-10CC	10cc	18cc	22cc	25cc	
	CS-20CC	20cc	36cc	44cc	50cc	

*Mat offers three bead sizes of 3mm, 4.5mm and 6mm

InterSep™ is manufactured by Pacific Bioceramics and distributed by Exactech, Inc.

Exactech, Inc. is proud to have offices and distributors around the globe.
For more information about Exactech products available in your country, please visit www.exac.com.