

Conquer TAA Logistics

VANTAGE[®] | Ankle
with **ACTIVE** Intelligence. | PSI

Patient-Specific Instrumentation Powered
by 3D Systems VSP[®] Solutions




exactech

Pre-plan resections, view deformities and determine implant sizing with the Vantage Ankle PSI tibia and talar cutting guides. Our 3D-printed instruments are designed for your patients' unique anatomies and tailored toward your specific surgical approach.

Accelerate to Surgery

Delivered within three weeks*

**Designed for a
Reproducible Fit**



Preserve the Periosteum

The instruments' unique soft tissue offsets are designed to allow you to preserve the periosteum while maintaining a reproducible fit.

Multiple cutting guides options to empower you to perform at your best.

Talar Cutting Guide




1. Flat Cut Talus Coupled Cutting Guide
2. Tibial Cutting Guide
3. Curved Talus Coupled Cutting Guide
4. Curved Talus Decoupled Cutting Guide
5. Flat Cut Talus Decoupled Cutting Guide

Cut Directly Through

Pre-plan tibia and talar resections to align and cut in one step.

**Tibia
Cutting
Guide**



To contact 3D Systems, please call 1-844-643-1001 to speak to a representative.

The Vantage Ankle PSI is manufactured by 3D Systems, Inc., and distributed only in the U.S. by Exactech.

*Orders will be fulfilled within three weeks as long as there are no extraordinary circumstances.