



EXACTECH|**BIOLOGICS**

Bone Marrow Aspirate Needle

Bone Marrow Aspirate Needle



FEATURES

- The needle includes five holes above the tip that allow for a more efficient aspiration.
- The stylet, with its pyramidal point, makes it possible to easily penetrate into the bone-marrow cavity.
- The ergonomically constructed handle is designed to enable a safer maneuverability, since the strength needed to penetrate the bone marrow cavity is distributed on the whole palm of your hand and not locally.
- The cannula is marked in 1cm increments that allow you to gauge the depth of insertion.

Code	Diameter/Length
RAN11	11G x 11cm
RAN1115	11G x 15cm
RAN11B*	11G x 11cm

* Supplied with an additional blunt stylet.

Exactech, Inc. is proud to have offices and distributors around the globe.
For more information about Exactech products available in your country, please visit www.exac.com.

Bone Marrow Aspirate Needles are manufactured by Ranfac Corp. and distributed by Exactech, Inc.