

PRECISION AND ACCURACY FOR **YOUR** • SHOULDER REPLACEMENT SURGERY



Using our Equinox Planning App before surgery, your surgeon can create a plan personalized to you. Paired with ExactechGPS during surgery, your plan will be performed with accuracy and precision. The Equinox Planning App and ExactechGPS are tools for your surgeon to use when making decisions for your surgery to allow for the best possible outcomes.



PLANNING

Equinox Planning App helps your surgeon create a personalized surgery plan.



PREDICTION

With Predict+, your surgeon can predict how you'll do after surgery.



EXECUTION

ExactechGPS assists your surgeon with precise placement of the implant during surgery.



PERSONALIZED PLANNING

The Equinox Planning App is an interactive computer program that creates a 3D model of your shoulder joint. This 3D model is developed from the CT scan you provided your surgeon. Your surgeon can select from implants considered for your surgery and connect them to the model of your shoulder joint. This allows your surgeon to make decisions such as what the best implant, implant size and placement is for you.



PREDICTION

Your surgeon will also have access to Predict+ through our planning app. Predict+ provides predictions for how you will do after surgery based on data collected about previous patients and surgeries. This data is analyzed using a type of artificial intelligence (AI) called machine learning, which finds patterns among previous patients and surgeries.

EXECUTION

Once your surgeon creates your plan in the planning app, it can be used with ExactechGPS - a visual computer navigation system that allows your surgeon to virtually navigate your joint and see exactly where to remove bone and implant your devices. ExactechGPS displays a 3D model of your shoulder joint and key measurements on a computer screen in the operating room. The visuals and measurements displayed on the screen are collected via a tracker placed on your shoulder during surgery. With real-time information and visuals about your shoulder joint and your joint movement, your surgeon can be accurate and precise with your implant placement.



Scan for more information.

From planning to the operating room, your surgeon is in full control. Our Equinox Planning App and ExactechGPS are software made to assist your surgeon with deciding what may be best for your surgery by providing data. However, your surgeon makes the final decisions about what will be best for you based on their expertise.

ExactechGPS, Equinox Planning App and Predict+ are developed by Blue Ortho, an Exactech subsidiary, and distributed by Exactech, Inc.